

## Congelatore professionale da Laboratorio -25°C



Ideale per la corretta conservazione a bassa temperatura.  
 Ottimo isolamento termico su tutte le pareti.  
 Interno in Acciaio Inox - Esterno verniciato  
 Ripiani estraibili e regolabili in altezza.  
 Porta con guarnizioni magnetiche e chiave.  
 Microinterruttore Apertura ante  
 Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C  
 Gestione degli Allarmi.  
 Sonde di temperatura al Platino PT1000.  
 Compressore ermetico con Gas Ecologico.  
 Alimentazione 230V 1f 50Hz

| CODICE               | GFREEZY Congelatore da Laboratorio                                                                                       | Dim Interne<br>LxPxH mm | Dim.Esterne<br>LxPxH mm | Peso<br>netto kg | Dim.Imballo<br>LxPxH mm | Peso<br>lordo kg |
|----------------------|--------------------------------------------------------------------------------------------------------------------------|-------------------------|-------------------------|------------------|-------------------------|------------------|
| <b>GFREEZ1400XC2</b> | GFREEZY 1400L Congelatore da Laboratorio -10-25°C, Ventilato, 4+4 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB | 1280x670x1500           | 1400x800x1990           | 229              | 1500x900x2170           | 249              |

## Professional Laboratory Freezer -25°C



Suitable for low temperature storage.  
Excellent thermal insulation on every wall.  
Stainless Steel internal chamber.  
Galvanized Steel exterior, PVC coated.  
Extractable and height adjustable shelves.  
Door with magnetic gaskets and lock with key.  
Microswitch at door opening.  
Control panel with digital temperature control system and 0,1°C display.  
Alarm control system.  
Platinum PT1000 temperature probes.

Hermetic compressor with ecofriendly gas.

Power supply 230V 1f 50Hz

| CODE          | GFREEZY Laboratory Freezer                                                                                                   | Interior Dim.<br>WxDxH mm | Exterior Dim.<br>WxDxH mm | Net<br>Weight kg | Packaging Dim.<br>WxDxH mm | Gross<br>Weight kg |
|---------------|------------------------------------------------------------------------------------------------------------------------------|---------------------------|---------------------------|------------------|----------------------------|--------------------|
| GFREEZ1400XC2 | GFREEZY 1400L Laboratory Freezer<br>-10-25°C, Ventilated, 4+4 Shelves,<br>Digital, Visual Acoustic Alarm,<br>DataLogger, USB | 1280x670x1500             | 1400x800x1990             | 229              | 1500x900x2170              | 249                |

