

## Doppio Frigorifero Professionale da Laboratorio



Ideale per la corretta conservazione a bassa temperatura.  
 Ottimo isolamento termico su tutte le pareti.  
 Interno in Acciaio Inox - Esterno verniciato.  
 Ripiani estraibili e regolabili in altezza.  
 Porta con guarnizioni magnetiche e chiave.  
 Porta vetrata per l'osservazione interna.  
 Microinterruttore Apertura ante.  
 Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C.  
 Gestione degli Allarmi.  
 Sonde di temperatura al Platino PT1000.

Compressore ermetico con Gas Ecologico.

Alimentazione 230V 1f 50Hz.

| CODICE        | GREFRY Doppio Frigorifero da Laboratorio                                                                          | Dim. Interne<br>LxPxH mm         | Dim.Esterne<br>LxPxH mm | Peso netto<br>kg | Dim.Imballo<br>LxPxH mm | Peso lordo<br>kg |
|---------------|-------------------------------------------------------------------------------------------------------------------|----------------------------------|-------------------------|------------------|-------------------------|------------------|
| GREFRY-300VFF | GREFRY Doppio Vano 2x150L<br>+2+12°C, Ventilato, 3+3 Ripiani,<br>Digitale, All.AcusticoVisivo,<br>DataLogger, USB | VS 480x465x720<br>VI 480x465x720 | 600x600x2040            | 125              | 700x700x2200            | 139              |

## Professional Laboratory Double Door Refrigerator



Suitable for low temperature storage.  
Excellent thermal insulation on every wall.  
Stainless Steel internal chamber.  
Galvanized Steel exterior, PVC coated.  
Extractable and height adjustable shelves.  
Door with magnetic gaskets and lock with key.  
Glass door for internal observation.  
Microswitch at door opening.  
Control panel with digital temperature control system and 0,1°C display.  
Alarm control system.  
Platinum PT1000 temperature probes.

Hermetic compressor with ecofriendly gas.

Power supply 230V 1f 50Hz

| CODE          | GREFRY Laboratory Double Door Refrigerator                                                                                  | Interior Dim.<br>WxDxH mm        | Exterior Dim.<br>WxDxH mm | Net<br>Weight kg | Packaging Dim.<br>WxDxH mm | Gross<br>Weight kg |
|---------------|-----------------------------------------------------------------------------------------------------------------------------|----------------------------------|---------------------------|------------------|----------------------------|--------------------|
| GREFRY-300VFF | GREFRY Double Compartment<br>2x150L +2+12°C, Ventilated, 3+3<br>Shelves, Digital, Visual Acoustic<br>Alarm, DataLogger, USB | UC 480x465x720<br>LC 480x465x720 | 600x600x2040              | 125              | 700x700x2200               | 139                |

